IKA

designed for scientists

IKA HBR 4 digital

IKA HBR 4 control

Betriebsanleitung	DE	04
Ursprungssprache		
Operating instructions	EN	15
Mode d'emploi	FR	26
Руководство пользователя	RU	37
使用说明	ZH	48
取扱説明書	JA	59
사용 설명서	KO	70

Indicaciones de seguridad	ES	81	Biztonsági utasítások	HU	91
Veiligheidsinstructies	NL	82	Varnostna opozorila	SL	92
Norme di sicurezza	IT	83	Bezpečnostné pokyny	SK	93
Säkerhetsanvisningar	SV	84	Ohutusjuhised	ET	94
Sikkerhedshenvisninger	DA	85	Drošības norādes	LV	95
Sikkerhetsanvisninger	NO	86	Nurodymai dėl saugumo	LT	96
Turvallisuusohjeet	FI	87	Инструкции за безопасност	BG	97
Normas de segurança	PT	88	Indicații de siguranță	RO	98
Wskazówki bezpieczeństwa	PL	89	Υποδείξεις ασφάλειας	EL	99
Bezpečnostní pokyny	CS	90			


Fig. 1

Pos.	Bezeichnung	Item	Designation
1	Istwert by Temperatur Istwert	1	Actual value by temperature actual value
2	Sollwert by Temperatur Sollwert	2	Target value by temperature target value
3	Phys. by physikalische	3	Phys. by physical
4	Betriebsart	4	Operating mode
5	Process value (Istwert)	5	Process value (actual value)
6	Set value (Sollwert)	6	Set value (target value)

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NAMUR commands		Function	Display (additional)
IN_NAME		Title request	
IN_PV_X	X=1;2;3;4;;	Current value reading	
IN_SOFTWARE		Software ID number, date, version request	
IN_SP_X	X=1;2;3;4;12;42 52;54;	Set target value reading	
IN_TYPE		Lab device identification request	
OUT_NAME name		Output of identification name. (Max. 6 characters; default: IKAHBR)	
OUT_SP_12@n		Setting WD safety limit temperature with set value echo	
OUT_SP_42@n		Setting WD safety limit speed with set value echo	
OUT_SP_X n	X=1;2;4;52; 54;	Setting of target value to n	
OUT_WD1@m		Watchdog mode 1: if event WD1 should occur, the heating and stirring functions are switched off and Er2 is displayed. Set watchdog time to m (20 - 1,500) seconds, with watchdog time echo. This command launches the watchdog function and must be transmitted within the set watchdog time.	
OUT_WD2@m		Watchdog mode 2: if event WD2 should occur, the speed target value is changed to the WD safety speed limit and the temperature target value is changed to the WD safety temperature limit value. The warning WD is displayed. The WD2 event can be reset with OUT_WD2@0 - this also stops the watchdog function. Set watchdog time to m (20 - 1,500) seconds, with watchdog time echo. This command launches the watchdog function and must be transmitted within the set watchdog time.	
RESET		Switches off the device function.	
RMP_CONT_X	X=1;4	Continuaton of ramp function (After prior RMP_PAUSE_X).	
RMP_IN_X	X=1;4	Reading the real segment number of ramp. With ramp not started: 0	
RMP_IN_X_y	X=1;4	Reading the accumulated value and the ramp segment duration of ramp segment y.	
RMP_LOOP_SET_X	X=1;4	Ending of ramp loop.	
RMP_LOOP_RESET_X	X=1;4	To work off the ramps in one loop.	
RMP_OUT_X_y n hh:mm:ss	X=1;4	Setting the accumulated value (n) and the ramp segment duration (hh:mm:ss) for ramp segment y.	
RMP_PAUSE_X	X=1;4	Stopping the ramp function. Freezing of real rated value and real ramp segment time.	
RMP_RESET_X		Switching off ramp functions and deleting of all set ramp segments	
RMP_START_X	X=1;4	Starting the ramp function, beginning with ramp function No.1 (Only possible after prior START_X. After RMP_STOP_X START_X is not necessary.)	
RMP_STOP_X	X=1;4	Switching off ramp function. Rated value = 0. (Ramp is maintained, that means, ramp can be restarted with RMP_STARTX_X).	
START_X</			

Maintenance and cleaning

The safety circuit must be checked over at least once a year by the user. To do this, fill the temperature control bath with 1 l of water and in operating mode A set the safety temperature to 100 °C. Then set the target temperature to 80 °C and start the heating function by pressing the Temp ON/OFF key. After the target temperature has been reached, switch the moderating bath off and back on again with the power switch and set the safety temperature to 70 °C. The temperature of the medium will then be 10 °C above the safety temperature, and the safety circuit will be activated. Error message Er22 or Er24 must be displayed.

The device is maintenance-free. It is subject only to the natural wear and tear of components and their statistical failure rate.

- Cleaning:
- For cleaning disconnect the mains plug!
 - Use only cleaning agents which have been approved by IKA to clean the devices:
 - Water containing surfactant / isopropyl alcohol.
 - Wear protective gloves during cleaning the devices.
 - Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.
 - Do not allow moisture to get into the device when cleaning.
 - Before using another than the recommended method for cleaning or decontamination, the user must ascertain with IKA that this method does not destroy the device.

- Spare parts order:
- When ordering spare parts, please give:
- device type
 - serial number, see type plate.
 - position number and description of spare part, see www.ika.com
 - software version.

Repair:

Please only send devices in for repair that have been cleaned and are free of materials which might present health hazards. For repair, please request the “**Safety Declaration (Decontamination Certificate)**” from IKA or use the download printout of it from IKA website at www.ika.com.

If your appliance requires repair, return it in its original packaging. Storage packaging is not sufficient when sending the device - also use appropriate transport packaging.

Accessories

For accessories see www.ika.com.

Warranty

According to **IKA’s** Terms and Conditions of sale and delivery, this product is covered by a warranty for a period of 24 months. In case of making a warranty claim, please contact your local dealer or, if you wish, you can send the device directly to our factory. Please include the sales invoice and state the reasons for your guarantee claim. In this case, you are responsible for shipping and handling costs.

The warranty does not cover wearing parts, nor to defects that are the result of improper use, insufficient care and maintenance or failure to follow the instructions in this operating manual.

Technical data

||
||
||